






ESTABLISHED 1907
WENIGER
EXCLUSIVE
MANUFACTURER OF

BUNNELL




WRIST-HAND-FINGER ORTHOSIS
STATIC-DYNAMIC – STATIC PROGRESSIVE POSITIONING
LATEX FREE-SOFT STRAPS-ORTHOSIS FITS RIGHT OR LEFT HAND



Expanded product line— Customer service
Please visit **BunnellSplints.com**
Sizing information on back page


BUNNELL ACTIVE HAND AND FINGER SPLINTS



#6

Spring Wire Safety Pin Splint to extend PIP Joint. Positional strap pad.


Sizes: XS, SM, MD, LG, XL



#6A

Modified Safety Pin Splint to extend PIP Joint. Dynamic Spring. Positional strap pad.


Sizes: XS, SM, MD, LG, XL



#6MINI

Mini Safety Pin Splint to extend DIP Joint. Positional strap pad.

Sizes: 1, 2, 3, 4



#6A MINI

Mini Modified Safety Pin Splint to extend DIP Joint. Dynamic Spring. Positional strap pad.

Sizes: 1, 2, 3, 4



#5

Finger Extension Splint w/Clockspring to extend a finger or thumb.


Sizes: XS, SM, MD, LG, XL




#5A

Finger Extension Splint w/Wrist Support. Dynamic Spring to M.P. joint for improved function. Malleable support extends to M.P.

Sizes: XS, SM, MD, LG, XL



#11
Finger Knuckle
Bender Splint to
flex the proximal
interphalangeal
joint.
Sizes: XS, SM,
MD, LG, XL



#11A
Reverse Finger
Knuckle Bender
Splint to extend the
proximal interpha-
langeal joint.
Sizes: XS, SM, MD,
LG, XL



#2
Knuckle Bender
Splint to flex the
metacarpophalangeal
joints.
Sizes: XS, SM, MD,
LG, XL




#2A
Knuckle Bender
Splint to flex the
metacarpophalangeal
joints w/outrigger to
extend the distal
joints.
Sizes: XS, SM, MD,
LG, XL




#2B
Knuckle Bender
Splint to flex the
metacarpophalangeal
joints with extra seg-
ment to flex the dis-
tal joints.
Sizes: XS, SM, MD,
LG, XL



#10
Composite Spring
and Elastic Splint to
extend the distal in-
terphalangeal joints.
Sizes: SM, MD, LG



#9
Reverse Knuckle
Bender Splint to
extend the metacar-
pophalangeal joints.
Sizes: XS, MD, LG,
XL



#9A
Reverse Knuckle
Bender Splint to
extend the metacar-
pophalangeal joints
with outrigger to
extend the distal
joints.
Sizes: XS, SM, MD,
LG, XL




#1
Spring Cock-up
Splint to dorsiflex
the wrist.
Kydex plastic
forearm shell.
Sizes: SM,
MD, LG



#1A
Spring Cock-up
Splint to dorsiflex
the wrist with
outrigger to flex
the M.P. joints.
Kydex plastic
forearm shell.
Sizes: SM, MD,
LG




#8
Dorsal Wrist
Splint to hold the
wrist in flexion
after repair of
flexor tendons or
after median and
ulnar nerve
suture. Kydex
plastic forearm
and wrist section.
Sizes: SM, MD,
LG



#8A
Dorsal Wrist
Splint to hold the
wrist in flexion
with outrigger to
extend distal
interphalangeal
joints. Kydex
plastic forearm
and wrist section.
Sizes: SM, MD,
LG




#4
Thomas Suspension
Splint to dorsiflex the wrist. For
radial palsy. Kydex plastic forearm shell. Outrigger
for thumb. Sizes: SM, MD, LG




#3
Modified Oppenheimer Splint Coiled Spring Wire to
dorsiflex the wrist for radial palsy. Thumb outrigger.
Hook and loop closure at wrist. Sizes: SM, MD, LG

#4.1
Thomas
Suspension
Splint for
cast or
thermo
plastic
application.
Sizes: One
size fits
most




#3.1
Modified Oppenheimer Splint with spring wire, to dorsiflex
the wrist for radial palsy. Spring wire applies less force to
extend wrist than #3. Thumb outrigger. Hook and loop
closure at wrist. Sizes: SM, MD, LG




#13

Combination of Coiled Spring Oppenheimer to dorsiflex the wrist with composite elastic splint to extend the distal joints. Thumb outrigger. Hook and loop closure at wrist. Sizes: SM, MD, LG




#13.1

Combination of Spring Oppenheimer to dorsiflex the wrist with composite elastic splint to extend the distal interphalangeal joints. Spring wire applies less force to extend wrist than #13. Thumb outrigger. Hook and loop closure at wrist. Sizes: SM, MD, LG



#13A

Combination of Coiled Spring Wire Oppenheimer to dorsiflex the wrist with Reverse Knuckle Bender to extend the metacarpophalangeal joints, and outrigger to extend distal joints. Thumb outrigger. Hook and loop closure at wrist. Sizes: SM, MD, LG




#13A.1

Combination of Spring Wire Oppenheimer with Reverse Knuckle Bender to extend the metacarpophalangeal joints and outrigger to extend distal joints. Spring wire applies less force to extend wrist than #13A. Thumb outrigger. Hook and loop closure at wrist. Sizes: SM, MD, LG



#13B

Combination of Coiled Spring Wire Oppenheimer to dorsiflex the wrist with Composite Elastic Splint to extend the M.P. and P.I.P. joints. Thumb outrigger. Hook and loop closure at wrist. Size: SM, MD, LG




#13B.1

Combination of Spring Wire Oppenheimer with composite elastic splint to extend M.P. and P.I.P. joints. Spring wire applies less force to extend wrist than #13B. Thumb outrigger. Hook and loop closure at wrist. Size: SM, MD, LG

STATIC PROGRESSIVE POSITIONING

- INCREASE R.O.M.—JOINT STIFFNESS—ELASTIC CONTRACTURES •
- COMFORTABLE PADDING • STAINLESS STEEL CONSTRUCTION •
- EASY TO USE HOOK & LOOP STRAPS TO CONTROL TENSION •



SIZE	LENGTH	
	INCHES	CM
XS	3 3/4	9.5
SM	3 7/8	9.8
MD	4 1/8	10.5
LG	4 3/8	11.1
XL	4 1/2	11.4

#18 P.I.P. Extension

Malleable finger platform and M.P. extension for improved function. Provided with optional use straps to secure distal portion of finger.



SIZE	LENGTH	
	INCHES	CM
XS	1 1/2	3.6
SM	1 5/8	4
MD	1 3/4	4.3
LG	1 7/8	4.8
XL	2	5.1

#19 P.I.P. and D.I.P. Flexion



#21 M.P. Extension with finger platform Shown with removable positional extension wedge attached. (Provided)

"Nonslip padding" provided for finger section.

Adjustable pads apply pressure to flex distal segment



Positional strap allows for optimal strap placement

#20 M.P. and P.I.P. Flexion with Distal Segment attached Hook and loop straps to allow positioning and removal of distal segment (above)


M.P. Flexion—Remove distal segment (not shown)

SIZING - #20, 21 & 28


SIZED ACCORDING TO WIDTH ACROSS M.P. JOINT

SIZE	SM	MD	LG
INCHES	2 1/2	3	3 1/2
CM	6.4	7.6	8.9

M.P. Extension straps



Wrist Extension Strap



#28 M.P. & WRIST EXTENSION

Separate straps to individually adjust MP and Wrist extension. Finger separators— "non-slip Padding" provided for finger section. Removable/positional wedge for increased MP extension included.

LOW PROFILE FINGER SPLINTS


- CLOSED CELL FOAM AND STAINLESS STEEL WIRE CONSTRUCTION
- EASY TO FIT AND ADJUST EXTENSION FORCE
- LATEX FREE
- CONTOURED PADS AND BANDS FOR IMPROVED FIT

SIZING— MEASURE LENGTH FROM DISTAL PALMER CREASE TO CREASE AT D.I.P.
POSITION PROXIMAL BAND OF SPLINT UNDER M.P. HEAD— DISTAL BAND BEHIND D.I.P.
JOINT

#24 PIP EXTENSION

Post injury—joint stiffness

	Length of splint	
FW01-1 XS	2 1/4"	5.7CM
FW01-2 SM	2 5/8"	6.7CM
FW01-3 MD	3"	7.6CM
FW01-4 LG	3 1/2"	8.6CM
FW01-5 XL	4"	10.0CM



.....FOR A LITTLE MORE FORCE

#25 Terminal PIP extension

	Length of splint	
FW07-1 XS	2 1/4"	5.7CM
FW07-2 SM	2 5/8"	6.7CM
FW07-3 MD	3"	7.6CM
FW07-4 LG	3 3/8"	8.6CM
FW07-5 XL	3 3/4"	9.5CM




#26 PIP SPRING EXTENSION—CAPNER TYPE DESIGN

Dynamic extension while allowing active flexion

Joint contractures that are elastic

	Length of splint	
FW02-2 SM	2 1/2"	6.4CM
FW02-3 MD	2 3/4"	7.0CM
FW02-4 LG	3"	7.6CM
FW02-5 XL	3 1/4"	8.3CM



CONTACT INFORMATION

951-654-2333 PHONE 951-654-2331 FAX
WWW.BUNNELLSPINTS.COM

SIZING INFORMATION

BUNNELL SPLINTS FIT RIGHT OR LEFT HAND

EASY TO ADJUST—STAINLESS STEEL CONSTRUCTION

•• *FINGER SPLINTS* ••

SIZED ACCORDING TO LENGTH OF SPLINT

MODEL # 5, 5A			MODEL # 6, 6A			MODEL # MINI 6, 6A			MODEL # 11, 11A		
SIZE	INCHES	CM	SIZE	INCHES	CM	SIZE	INCHES	CM	SIZE	INCHES	CM
XS	6	15.2	XS	2	5.1	1	1 5/16	3.3	XS	1 1/2	3.6
SM	6 1/2	16.5	SM	2 1/2	6.4	2	1 9/16	4	SM	1 5/8	4
MD	7	17.8	MD	3	7.6	3	1 7/8	4.8	MD	1 3/4	4.3
LG	7 1/2	19.1	LG	3 1/2	8.9	4	2	5.1	LG	1 7/8	4.8
XL	8	20.3	XL	4	10.2				XL	2	5.1

•• *HAND FINGER & WRIST HAND FINGER SPLINTS* ••

SIZED ACCORDING TO WIDTH ACROSS M.P. JOINTS


HAND FINGER SPLINTS

MODEL # 2, 2A, 2B, 9, 9A, 10, 10A

WRIST HAND FINGER SPLINTS

MODEL # 1, 1A, 13, 13.1, 13A, 13A.1,
13B, 13B.1

SIZE	WIDTH ACROSS M.P. JOINT
XS	2"
SM	2 1/2"
MD	3"
LG	3 1/2"
XL	4"



WENIGER

915 INDUSTRIAL WAY

SAN JACINTO, CA 92582-3890

PH. (951) 654-2333 FAX. (951) 654-2331

WWW.BUNNELLSPINTS.COM